Non-contact Ultrasonic Thermal Imaging Module

Introduction

Adaptive multi-sensor infrared temperature measurement algorithm.

The Smart temperature terminal achieves accurate infrared temperature measurement with data display of body temperature, voice alarm. Fast and accurate body temp measurement with long-distance detection to minimize close interactions. This system runs completely autonomous to allow temperature screening with less human interaction.

LTK05T



- Built-in non-contact infrared temp measurement module to ensure accurate temp measurement.
- Temperature measurement range is from 35 °C to 42 °C with accuracy ± 0.3 °C.
- · Self-service temperature measurement to avoid cross infection
- Alarm notification when the temperature of the person being detected is abnormal.
- Advanced aluminum alloy metal housing and durable.
- Support wall and pole installation

Applications



Commercial building/ Community



Cinema



School/ Camp



Bank/Hotel

Specifications


| Item | Specification | |
|-------------------------|--|--|
| CPU | ARM Cortex-m | |
| Temperature Display | LED digital display | |
| Sensor | Thermal image infrared temperature sensor | |
| Temperature measurement | Range: 35°C ~42°C , accuracy: ±0.3°C | |
| Detection distance | 30-60cm | |
| Respond time | 1~3s | |
| Alarm method | Warning sound (70db)+ red beating temp indicator | |

Specifications

| General specifications | | |
|------------------------|----------------------------|--|
| Power consumption | ≤5W | |
| Working temperature | 15°C ~35°C (indoor) | |
| Working humidity | 0%-85% RH (Non-condensing) | |
| Dimensions (L×W×H) | 45mm×105mm×193mm | |
| Weight | 542g | |

* Specifications are subject to change without notice.

Descriptions



Installation

For wall

- 1 About ten seconds after the power is turned on, the device starts to complete;
- 2 To fix the device to wall through the back mounting point.
- 3 To ensure the device screen in front of human face when it's installed above or to the side of the passage.
- 4 The recommended angle to install is 0-20degrees.
- 5 To read body temp measurement correctly, the recommended distance between the device and collection point is 0 30-60 cm.





For pole

- 1 About ten seconds after the power is turned on, the device starts to complete;
- $2\,-\,$ To fix the device to pole through the bottom mounting point.
- 3 Set up the device and pole directly in front of the passage with display facing the human face.
- 4 The recommended installation height is 140-160cm.
- 5- To read body temp measurement correctly, the recommended distance between the device and collection point is 30-60 cm.







- Device should be installed indoors, make sure there is no wind between device and person.
- No direct sunlight to avoid incorrect temperature measurement caused by wind dissipation and direct sunlight.

Packing List

| Item | Model | QTY | Note |
|---------------|----------|-----|---|
| Machine body | JRI . 3R | 1 | |
| Bracket | 1 | 1 | |
| Adapter | 1 | 1 | 5V1A AC-DC power adapter |
| Charging line | 1 | 1 | 1 meter, USB Type-A to USB Type-C charge wire |
| User manual | 1 | 1 | |